



Zepto Life Technology, Inc.

MEDICAL TEST MENU

Effective Date: 17 July 2026

FUNGAL

Test Code	Test Name	CPT Code	Specimen Requirements		Shipping	Rejection Criteria	TAT
FM001	FungiFlex® Mold Panel	87798	Acceptable Specimen: K2-EDTA whole blood or plasma	Whole Blood Specimen Stability: Whole blood collected in an K2-EDTA tube may be stored at room temperature for up to 24 hours before plasma separation is required.	Whole Blood Shipping: Ship at ambient temperature only if Zepto will receive the specimen within 24 hours of collection . If delivery will exceed 24 hours, separate plasma before shipment. Plasma Shipping: If Zepto is expecting to receive the specimen within 96 hours of collection , ship at ambient temperature. If receipt is expected beyond 96 hours , freeze the plasma ($\leq -20^{\circ}\text{C}$) and ship on dry ice .	Specimens missing two unique patient identifiers will be rejected.	24 hours*
FPJ01	FungiFlex® <i>Pneumocystis jirovecii</i>		Whole Blood Preferred volume: ≥ 6 mL Minimum volume: 4 mL Plasma Preferred volume: ≥ 3 mL Minimum volume: 2 mL			Leaking, broken, or damaged specimen container. Gross hemolysis, clotting, significant cellular contamination, or specimen deterioration. Specimens collected in an expired collection tube. Specimens shipped under unacceptable conditions (e.g., unfrozen when freezing was required).	
FMPJ	FungiFlex® Mold Plus Panel		Acceptable Specimen: K2-EDTA whole blood or plasma Whole Blood Minimum volume: 8 mL Plasma Minimum volume: 4 mL	Plasma may be stored frozen at $\leq -20^{\circ}\text{C}$ for up to 4 weeks. Note: Plasma specimens should not undergo more than three freeze-thaw cycles.		Note: Improper specimen collection, inadequate plasma separation, hemolysis, clotting, cellular contamination, specimen deterioration, contamination, or improper storage and transport conditions may reduce target DNA recovery, affect assay performance, or result in specimen rejection.	

Additional Information:

*Results are reported within 24 hours of specimen receipt at Zepto Monday - Friday. Requested retesting and holidays may extend TAT.

Disclaimer: FungiFlex® test offerings are laboratory-developed tests (LDT) that have not been cleared or approved by the U.S. Food and Drug Administration (FDA). These tests have been validated by Zepto Life Technology in accordance with CMS CLIA regulations to establish its performance characteristics.